

IPS 14


infusion mode

IPS
SERIES

// 14 IPS ➤

MODEL	IPS14
MODE	Infusion only
SYRINGES	
Capacity	2 Independent
Type	Glass, plastic, steel
Min Size	Hamilton 0.5 µl (ID:0.103mm)
Max Size	Monoject 140 ml (ID:38mm)
FLOW RATE	
Min	10 nl/min (Hamilton 0.5 µl)
Max	120 ml/min (140 ml Monoject Syringe)
MOTOR TYPE	Stepper Motor
MOTOR DRIVE CONTROLLER	Microcontroller with Microstep Drive
STEP RESOLUTION (Travel/Microstep)	357 nm/µstep
STEP RATE	
Min	10 sec/µstep
Max	200 µsec/µstep
LINEAR SPEED RATE (PUSHER TRAVEL RATE)	
Min	2,14 µm/min
Max	107 mm/min
LINEAR FORCE	20 kg for each syringe
CONNECTORS	Wi-fi (optional)
RECIPE SAVE-RECALL	No
POWER SUPPLY	12 VDC 1A (110-240 VAC adapter included)
DIMENSIONS	W:200 L:257 H:137 mm
WEIGHT	3.3 kg
ENVIRONMENTAL SPECIFICATIONS	
Operating Temperature	Between -10 and 50 (°C)
Storage Temperature	Between -30 and 85 (°C)
Operating Humidity	RH 10% - RH 90% (not condensed)
Storage Humidity	RH 10% - RH 90% (not condensed)

- 2 Independent channels
- 4.3" Resistive Touch TFT LCD screen
- Microcontroller software
- User-friendly interface
- High precise pumping
- Adjustable flow rate even during process



IPS